

# A New Convenient System



## **Disposable Proctorectoscope Tubes**

Packed sterile – for immediate use

# Disposable Proctoscope Tubes

## A New Convenient System

Richard Wolf Disposable Proctoscope Tubes have a long track record of providing a clear view in proctological diagnostics for many years.

This step has now become even more efficient. Sterile packing at the factory means that the proctoscope tubes are now ready for use immediately. This eliminates the process of disinfection prior to application.



### ■ Improved handling

Each tube is packed singly and sterile. Disinfection prior to application is eliminated!

### ■ Cost-effective solution

You benefit from the ideal combination of disposable tube and reusable handle.

### ■ Everything in view

An easily readable scale makes documentation and checks easier.

### ■ Everything under control

The reusable handle permits easy application with optimum illumination.

#### Tubes and accessories

	<b>Rectoscope tube</b> , sterile, Ø 21 mm, working length 255 mm, pack of 10 .....4844255
	<b>Proctoscope tube</b> , sterile, Ø 21 mm, working length 120 mm, pack of 10 .....4844120
	<b>Handle</b> with annular cold light illumination .....8838.90
	<b>Hand bulb, sterile (pack of 5)</b> Hand bulb with filter for disposable use Ø 49 mm, OL 567 mm .....412700
	<b>Locking cap</b> with observation window .....8838.60

